SOEHNLE

DE

EN

FR

IT

NL

ES

SV

CS

PL





SYSTO MONITOR CONNECT 400

GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MODE D'EMPLOI
ISTRUZIONI PER L'USO
GEBRUIKSAANWIJZING
INSTRUCCIONES DE MANEJO
BRUKSANVISNING
NÁVOD K POUŽITÍ
INSTRUKCJA OBSŁUGI

 Bluetooth®

Deutsch	2
English	18
Français	32
Italiano	48
Nederlands	64
Español	80
Svenska	96
Český	112
Polsku	128

EINFÜHRUNG

DE

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Guangdong Transtek Medical Electronics Co., Ltd.
Zone A, No.105, Dongli Road, Torch Development District, Zhongshan, 528437, Guangdong, China

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MDSS - Medical Device Safety Service GmbH,
Schiffgraben 41,
30175 Hannover/Germany

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
INTRODUCTION

EN

INTRODUCTION.....	18	DATA MANAGEMENT	26
Safety information		Recalling stored values	
Scope of delivery		Deleting stored values	
General description		Data transmission	
Indications for Use		INTERPRETING THE RESULTS.....	27
Measurement Principle		CUSTOMER INFORMATION	28
Device description		Maintenance and Storage	
Power supply		Contact	
Inserting and replacing the batteries		Customer Service	
The LCD screen		Guarantee	
PRIOR TO THE FIRST MEASUREMENT..	24	TROUBLESHOOTING	29
Setting date, time, and unit of measure		TECHNICAL DATA.....	30
Selecting the user		LIST OF COMPLIED EUROPEAN STANDARDS	31
Soehnle Connect app		EMC GUIDANCE	144
MEASURING.....	25		
Measurement information			
Cuff placement			
Starting the measurement			

Safety information

	Symbol for "THE OPERATION GUIDE MUST BE READ"		Symbol für „MIT HOCHFREQUENZSENDER“
	Symbol for "COMPLIES WITH MDD 93/42/EEC REQUIREMENTS"		Symbol for "TYPE BF APPLIED PARTS"
	Symbol for "MANUFACTURER"		Symbol for "ENVIRONMENT PROTECTION - Electrical waste products should not be disposed of with household waste. Please recycle where facilities exist. Check with your local authority or retailer for recycling advice"
	Symbol for "SERIAL NUMBER"		Symbol for "Authorised Representative in the European Community"
	Symbol for "DIRECT CURRENT"		Caution: These notes must be observed to prevent any damage to the device.
	The Bluetooth Combination Mark		

**IMPORTANT INFORMATION! PLEASE KEEP!**

Carefully read these user instructions and in particular the safety information prior to using this device and keep these user instructions for future use. When handing the device over to third parties, please also include these user instructions. Non-compliance with these instructions may cause injuries or damage the device.

**Information on use**

- The device is intended for the independent measurement and monitoring of arterial blood pressure. Only a medical doctor can provide a diagnosis.
- The user must review the functionality and proper condition of the device prior to using it. Continued use of a damaged device may result in injuries and erroneous measurement results or pose a serious risk. In case of technical problems with the device, please contact Soehnle's CUSTOMER SERVICE. Never open or repair the device.
- The device is intended to be used only by adults. Do not use the device if you are pregnant or suspect a pregnancy.
- This device may be used only for the purpose described in these instructions. The manufacturer does not accept any liability for damages based on improper use.
- The device is designed only for personal use.
- It must not be used on any extremities but the arms.
- If you take medication please ask your doctor about the best time to measure your blood pressure.
- If the device is used for measurements on patients with widespread arrhythmia, e.g. extrasystoles of the atrium or ventricularextrasystoles, the best results are generally achieved with deviations. Please discuss the results with your doctor.
- If the cuff pressure is higher than 40 kPa (300 mmHg), the device releases air automatically. If the cuff does not release air at a pressure of 40 kPa (300 mmHg), please remove the cuff from your arm and push the START/STOP key to prevent further inflation.
- The device is not approved for AP or APG categories.
- Avoid interference signals emitted by strong electromagnetic fields or rapid transient disturbances to prevent measuring errors.
- The device is not suited for continuous blood pressure monitoring in medical emergencies. Blood circulation may be disrupted.
- The cuff material has been tested and is proven to correspond to the requirements of ISO 10993-5:2009 and ISO 10993-10:2010. It does not cause skin irritation or other reactions.
- Please use only the ACCESSORIES and removable parts prescribed

INTRODUCTION

EN

by the manufacturer. Otherwise, danger may arise for the user and the device.

- Please inform Soehnle customer service of unexpected operating status or events.
- The optional adapter is a component of the MEDICAL ELECTRICAL DEVICE.
- The plugs/pins of the adapter plug isolate the device against the mains power supply.
- Please do not place the device in a location in which it is difficult to disconnect from the power supply.
- Keep the device outside of the reach of small children/pets to prevent the inhalation or swallowing of small parts. The cable/cord poses a risk of strangulation.
- When stored at the minimum storage temperature between uses, electrical medical devices must be warmed up for at least 30 minutes before they are operational. When stored at the maximum permissible storage temperature between uses, electrical medical devices must cool off for at least 30 minutes before they are operational.
- Please carefully monitor a possible interruption of blood circulation and a possible injury of the patient.
- Do not place the cuff over a wound as this may cause further injuries.
- Do not place the cuff on the same arm on which other electrical medical devices are placed.

- Please ensure that the operation of the device does not interrupt blood circulation for a longer period.
- Avoid squeezing or constriction of the connecting tube.
- The Blood Pressure Monitor has been clinically tested pursuant to the requirements of EN 1060-4:2004.

Metrological maintenance and care information

- Do not perform any metrological maintenance if the device is in use.
- For reliable operation, we recommend that the device is subject to metrological control once every two years pursuant to Sec. 14 of the Medical Product Operating Directive (MPBetriebV). We recommend testing the device's performance once every two years and after metrological control and repair, at a minimum by testing the tolerance requirements of the pressure display and tightness of the cuff (testing at least at 50 mmHg and 200 mmHg).
- The metrological control must be performed by a licensed tester or testing institute.
- We recommend cleaning the device and cuff only with a soft cloth.

Storage information

- If the device is not being used, store it and the adapter in a dry location away from heat, dust, fibres, lint, and out of direct sunlight. Do not place heavy objects on the storage box.

 **Information about handling batteries**

- Do not mix new and used batteries.
- Do not use different types of batteries at the same time.
- Do not throw batteries into a fire. Batteries may explode or leak. Remove the batteries if the device is not used for a long time.

 **Disposal information**

**Battery disposal
EC Directive 2008/12/EC**

Batteries are not a part of your regular household waste. You must return batteries or rechargeable batteries to your municipality's public collection or wherever batteries of the respective type are being sold.

**Disposal of electric and
electronic devices, EC
Directive 2012/19/EU**

This product is not to be treated as regular household waste but must be returned to a collection point for recycling electric and electronic devices. Further information is available from your municipality, your municipality's waste disposal services, or the retailer where you purchased your product.

CE **Declaration of conformity**
0123 This device is in compliance with the basic requirements and the other pertinent regulations of the Directive 93/42/EEC.

EN

Scope of delivery

- 1x blood pressure monitor Systo Monitor Connect 400
- 1x cuff (22 – 42 cm / 8.66 – 16.54 inch)
- 4x alkaline battery AAA 6 V DC
- 1x user instruction

A power adapter is NOT included in the scope of delivery. It is available separately, if needed. Information is available under "Power supply".

General description

Thank you for choosing the blood pressure monitor Soehnle Systo Monitor Connect 400 for use on the upper arm.

Measurements taken with Soehnle Systo Monitor Connect 400 are equal to those taken by a trained professional.

This user instruction contain important information concerning safety, care of the device, and explain the proper use step-by-step.

INTRODUCTION

Indications for Use

EN The Soehnle Systo Monitor Connect 400 is digital monitors intended for use in measuring blood pressure and heartbeat rate with arm circumference ranging from 22 – 42 cm / 8.66 – 16.54 inch.

It is intended for adult indoor use only.

Measurement Principle

This product uses the Oscillometric Measuring method to detect blood pressure. Before every measurement, the unit establishes a “zero pressure” equivalent to the air pressure. Then it starts inflating the arm cuff, meanwhile, the unit detects pressure oscillations generated by beat-to-beat pulsatile, which is used to determine the systolic and diastolic pressure, and also pulse rate.

Device description

Please note the fold-out page.

- A** Cuff
- B** Air hose
- C** LCD display
- D** Pump
- E** Power plug connection
- F** Battery compartment
- G** Memory key (Q)
- H** Settings key (⊗)
- I** Start/Stop key (⏻_{START})

Power supply

1. Battery operation:

4 alkaline batteries AAA 6 V DC

2. Operation with optional power adapter:

6V \equiv 1A

Please use only the adapter approved by the manufacturer. Available separately (Item No. 89165/5).

Inserting and replacing batteries

1. Open the battery compartment.
2. Insert the batteries. Note the correct polarity.
3. Close the battery compartment cover.

The normal service life of new and unused batteries is 300 measurements at a measurement period of 60 seconds each.

The LCD screen

Symbol	Description	Explanation
SYS	Systolic blood pressure	Upper blood pressure value
DIA	Diastolic blood pressure	Lower blood pressure value
PUL /min	Pulse display	Pulse in beats per minute
	Symbol pressure release	Air is released from the cuff
	Memory	Displays that the device is in memory mode
kPa	kPa	Unit of measurement for blood pressure (1 kPa = 7.5 mmHg, see page 24)
mmHg	mmHg	Unit of measurement for blood pressure (1 mmHg = 0.133 kPa, see page 24)
	Battery charge display	Displays charge status of the battery
	Irregular pulse	Irregular pulse detected during measurement
	Blood pressure level	Indicates blood pressure level (based on WHO scale, see page 27)
	Current time	Year/month/day, hour/minute
	Heartbeat	The blood pressure monitor detects the heartbeat during the measurement
	User 1	Active user profile (user 1)
	User 2	Active user profile (user 2)
	Movement indicator	Indication of result distorting movement
	Average value	Average value of the blood pressure
	Bluetooth symbol	Shows that Bluetooth is activated
	Data transmission	Shows that no data was transmitted to the smartphone

PRIOR TO THE FIRST MEASUREMENT

EN

Setting date, time, and unit of measure

Date and time must be set prior to using the blood pressure monitor so that every measurement value that is stored can be time stamped.

1. If the device is switched off, push the key (⊙) to display the time. Then push and hold the key (⊙) to set the year.
2. Push the key (Q) to change the [YEAR]. Each time you push the key, the year increases by one.
3. Once the year has been set, push key (⊙) to save the information and continue.
4. Repeat steps 2 and 3 to set [MONTH], [DAY], [HOUR], [MINUTE] and [UNIT OF MEASURE]

You can switch between the unit of measure "kPA" (kilo pascal) or "mmHg" (millimetre mercury column), where mmHg is the most common unit for measuring blood pressure.

5. After setting the unit of measure, the LCD display shows "done" and then displays all settings entered by you before switching off.

Selecting the user

A total of two users can be programmed. 60 datasets can be stored for each user.

1. If the device is switched off, push and hold the key (Q) to select the user. The user symbol (♂) flashes.
2. Push key (Q) again to select between user 1, 2 or "Guest" (guest mode).

The measurements are not stored or transmitted in guest mode.

3. Confirm the selection of the respective user by pushing (⊙). The respective user and "done" are shown before the display switches off.

Soehnle Connect app

To use all functions of your Systo Monitor Connect 400, we recommend that you connect the device to the Soehnle Connect app prior to first use.

1. Download and install the free Soehnle Connect app on your smartphone or tablet. Scan the QR code to access the app directly. Or find it in Google Play or the App Store..

Android: min. 4.4.4 / iOS: min. 9.0 / Bluetooth: min. 4.0

2. Activate Bluetooth on your smartphone or tablet and open the Soehnle Connect app. When opening the app for the first time, you will automatically be guided through the process to connect a Soehnle end device.

Once you have already used the Soehnle Connect app, to connect in the app go to the "main menu" and press "my devices" and then "add devices".

Compatibility

The Soehnle Connect app is compatible with Google Fit® and Apple Health®.

Measuring information

The measurements may be incorrect under the following circumstances:

- Within 1 hour after a meal or a drink
- Measurements immediately after consuming tea, coffee, or after smoking
- Within 20 minutes after a bath
- While speaking or moving fingers
- In a very cold environment
- With a strong urge to urinate
- Relax for 5 min. before measuring.
- Wait at least 3 min. between two measurements. This allows the blood flow in your arm to return to normal.
- For a reasonable comparison of various measurement values, always attempt to measure in the same situations. For example, take your measurement daily at around the same time and with the upper arm in the position as shown to you by a doctor.

Cuff placement

1. Place the cuff on the upper arm and then turn the hose towards the inside of the arm in line with the pinkie finger. You can also align the arterial mark so that it is placed above the largest artery (on the inside of the arm). You can find the largest artery by gently pushing 2 fingers approximately 2 cm above the crook of your arm, on the inside of your left arm. Find where the pulse can be felt the strongest. This is where the largest artery is located. **(Fig. 1)**

2. The cuff should be snug but not too tight. It must be possible to insert one finger between the cuff and your arm. **(Fig. 2)**
3. Sit in a comfortable position and place the arm on which the measurement will be taken on a level surface.

Sit down comfortably prior to starting the measurement: do not cross your legs, place the soles of your feet on the floor, your back and arm should be supported.

4. For users with hypertension:

The centre of the cuff should be at the level of the right atrium **(Fig. 3)**.

Starting the measurement

1. If the blood pressure monitor is switched off, switch it on by pushing the key (). The device performs the measurement.
 - The device will first zero adjust itself
 - The cuff will then be inflated and the device will measure the blood pressure and pulse.

Push the key () to stop the measurement at any time.

- Once the measurement is completed, the air is released from the cuff and the measured values will be displayed. They will be automatically allocated to the selected user and stored. Year, date, and time are shown alternately.
2. Push the key () to switch the device off. The device switches off automatically after 1 minute.

Recalling stored values

EN

1. Push (Q) to switch the blood pressure monitor on. The average value of the last three datasets will be displayed.
2. Push the key (Q) or (O) to get to the desired datasets. While scrolling, the number of the displayed dataset and the total number of stored measurements will be shown in the upper right corner.

If you leave a dataset displayed, date and time of the measurement will also be shown. They will be put through automatically.

3. If you would like to recall the datasets of another user, first push the key ($\frac{\text{O}}{\text{START}}$) to switch the device off if the blood pressure monitor is still in the “recall memory” mode.

Push and hold the key (Q) to enter the mode “select user”, push the key (Q) again to select between user 1 and user 2. Push the key (O) to confirm the user, then recall the measured datasets of the selected user with (Q).

The most current dataset (1) is shown first. Each new measurement is stored at dataset (1). All other datasets are moved back by one (2 becomes 3, etc.) and the last dataset (60) is removed from the list.

Deleting stored values

If a measurement is not performed correctly, you can delete all measurement results for the selected user as follows.

1. Push and hold key (Q) for approximately 3 seconds while the blood pressure monitor is in the “recall memory” mode. “dEL ALL” (delete ALL) will flash and the selected user will be displayed.

Push ($\frac{\text{O}}{\text{START}}$) to exit the “delete” mode without deleting datasets.

2. Push (O) to confirm deletion, and “user+dEL done” will be displayed. The blood pressure monitor then switches off automatically.

Data transmission

After each measurement, your personal blood pressure information is transmitted via Bluetooth to your mobile device on which the Soehnle Connect app is installed.

1. After the measurement, the Bluetooth symbol (B) flashes on the top right of the screen. If a connection with your mobile device is established, the Bluetooth symbol (B) freezes.
2. The new measurement data will be transmitted after a successful connection with your mobile device has been established.
3. After successful transmission, the Systo Monitor Connect 400 displays “done” on the screen.

In the event of an error during data transmission, all stored values (for both users) can be re-send. To do this, push and hold the start/stop button for 3 seconds until the Bluetooth symbol flashes on the screen. Then repeat steps 2 and 3.

What is systolic and diastolic blood pressure?

When the heart ventricles contract to pump blood out of the heart, the blood pressure reaches its maximum value in the cycle, the so-called systolic value. The blood pressure reaches its lowest value in the cycle when the heart ventricles relax. This is the diastolic blood pressure.

How are blood pressure values usually categorized?

The classification of blood pressure published in 1999 by the World Health Organization (WHO) and the International Society of Hypertension (ISH) is as follows:

Classification	SYS (mm HG)	DIA (mm HG)
Grade 3 HT	≥ 180	≥ 110
Grade 2 HT	160 – 179	100 – 109
Grade 1 HT	140 – 159	90 – 99
High normal	130 – 139	85 – 89
Normal	120 – 129	80 – 84
Optimal	< 120	< 80

HT = Hypertension (high blood pressure)

IMPORTANT

Only a medical doctor can determine your normal blood pressure range. Please see a doctor if your measurement results are outside of this range.

Recognizing irregular pulse

The device recognizes irregular pulse (arrhythmia) if the heartbeat fluctuates while measuring the systolic and diastolic blood pressure. The blood pressure monitor records heartbeat intervals and calculates the average value during each measurement. If an average value is greater than or equal 25%, the measurement result also displays a warning signal for irregular pulse.

The warning “irregular pulse” shows that an irregular pulse has been determined during the measurement, which might indicate arrhythmia. Generally, this is NOT a cause for concern. However, if the symbol is displayed frequently, we recommend that you speak with your doctor. The device does not replace a cardiology examination but provides an early detection of irregular pulse.

Maintenance and Storage

EN

Your device will provide the best measurement accuracy if you comply with the following instructions:

- Store in a dry location and avoid direct sunlight
- Avoid contact with water, clean with a soft cloth if necessary.
- Avoid strong vibrations and impact.
- Do not store in dusty environments or varying temperatures.
- Do not clean the reusable cuff with water and do not submerge it in water.

Contact

Soehnle Customer Service

Hotline +49 (0) 26 04 – 977 0
Mon – Fri,
7:00 a.m. – 6:00 p.m.

Further information about our products can be found at www.soehnle.com.

Distributor:

Soehnle is a brand of Leifheit AG,
Leifheitstraße 1
56377 Nassau/Germany

Manufactured by:

Guangdong Transtek Medical
Electronics Co., Ltd.
Zone A, No.105, Dongli Road, Torch
Development District, Zhongshan,
528437, Guangdong, China

European representative:

MDSS - Medical Device Safety
Service GmbH, Schiffgraben 41,
30175 Hanover, Germany

Guarantee

Leifheit AG offers a 2 year guarantee for this quality product as of the purchase date (or in case of orders, as of receipt of the goods). Guarantee claims must be asserted within the guarantee period and immediately after occurrence of the defect. The guarantee includes product quality.

The following are excluded from the guarantee:

1. Wear caused by use or other natural wear,
 2. Damage caused by improper use or handling (e.g. impact, shock, drop),
 3. Damage caused by non-compliance with these user instructions.
- In case of a guarantee claim, Leifheit, at its own discretion, either repairs the defective components or replaces the product. If a repair cannot be performed and an identical product is no longer available in the product range for the purpose of replacement, you will receive a replacement product of the closest equivalent value. Reimbursement of the purchase price is not possible in case of a guarantee claim. This guarantee furthermore does not grant any claims for damages. To assert the guarantee, please contact the seller where you have acquired the product and present the defective product and the purchase receipt (copy). This guarantee is valid worldwide. Your statutory rights, in particular warranty rights, continue to apply and are not limited by this guarantee.

PROBLEM	POSSIBLE CAUSE	RESOLUTION
Display does not light up	Batteries are empty	Replace batteries
	Batteries are not inserted correctly.	Insert batteries correctly.
	The DC adapter is not inserted correctly.	Insert the DC adapter correctly.
Lo +  is displayed	Low battery charge	The charge status of the battery is too low; batteries must be replaced.
Display is difficult to see or slow	Battery charge is too low.	Replace batteries
E 1 is displayed	The cuff is too tight or too loose.	Reposition cuff and repeat measurement
E 3 is displayed	The pressure of the cuff is excess.	Relax for a moment and repeat the measurement.
E 10 or E 11 is displayed	The blood pressure monitor displays movement during the measurement.	Movements may affect the measurement. Relax for a moment and repeat the measurement.
E 20 is displayed	No pulse was detected during measurement.	Remove the clothing from your arm and repeat the measurement.
E 21 is displayed	A precise measurement could not be performed.	Relax for a moment and repeat the measurement.
EExx is displayed.	A "calibration error" occurred. (XX can indicate a digital signal, e.g. 01, 02 etc.)	Repeat the measurement. If the problem occurs again please contact your seller or our customer service for further assistance.
"Out" is displayed	Outside of the measurement range	Relax for a moment. Reposition the cuff and repeat the measurement. If the problem occurs again please contact your doctor.

TECHNICAL DATA

Item number	68097	(Modell: TMB-1585-BT)
Power supply	Battery operation: 4 alkaline batteries AAA 6V DC Power adapter operation 6V \equiv 1 A (Please use only the recommended power adapter.)	
Display mode	LCD display (100 x 68 mm / 3.94 x 2.68 inch)	
Measurement mode	Oscillographic testing mode	
Measurement range	Target cuff pressure: 0 mmHg – 300 mmHg (0 kPa – 40 kPa) measurement pressure SYS: 60 mmHg – 230 mmHg (8.0 kPa – 30.7 kPa) DIA: 40 mmHg – 130 mmHg (5.3 kPa – 17.3 kPa) Pulse: (40 – 199) beats per minute	
Accuracy	Pressure: 5 °C – 40 °C within \pm 3 mmHg (0.4 kPa) pulse at: \pm 5 %	
Normal operating conditions	Temperature: 5 °C to 40 °C Relative humidity: \leq 85 % RH Atmospheric pressure: 86 kPa to 106 kPa	
Storage and transport conditions	Temperature: -20 °C to 60 °C Relative Humidity: 10 % RH-93 % RH Atmospheric Pressure: 50 kPa-106 kPa	
Measurement perimeter on the upper arm	22 – 42 cm / 8.66 – 16.54 inch	
Weight	approx. 277 g / 9.7 oz (excluding batteries and cuff)	
External dimensions	approx. 154 x 106 x 57 mm / 6.06 x 4.17 x 2.24 inch	
Accessories	4 alkaline batteries AAA	
Operating mode	Continuous operation	
Protection category	Applied part type BF	
Water ingress protection	IP21 means that the blood pressure monitor is protected against solid foreign bodies with a size of 12.5 mm / 0.49 inch and larger as well as vertically dripping water.	
Device Classification	Battery Powered Mode: Internally Powered ME Equipment AC Adaptor Powerde Mode: Class II ME Equipment	
Software version	V01	
Bluetooth	Bluetooth Module No.: AW 2540 MV 1 High frequency range: 2402 MHz to 2480 MHz Output power level: -1 dBm Supply Voltage: 2V - 3.6V Transmitting Distance: 10 meters	

WARNING: Modifications of the device are not permitted.

LIST OF COMPLIED EUROPEAN STANDARDS

EN

Risk management	EN ISO 14971:2012 Medical devices - Application of risk management to medical devices
Labeling	EN 980:2008 Symbols for use in the labelling of medical devices
User manual	EN 1041:2008 Information supplied by the manufacturer of medical devices
General Requirements for Safety	EN 60601-1:2006 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance EN 60601-1-11:2015 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment
Electromagnetic compatibility	EN 60601-1-2:2007 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
Performance requirements	EN ISO 81060-1:2012 Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type EN 1060-3:1997+A2:2009 Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems
Clinical investigation	EN 1060-4:2004 Non-invasive sphygmomanometers - Part 4: Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers
Usability	EN 60601-1-6:2010 Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability EN 62366:2008 Medical devices - Application of usability engineering to medical devices
Software life-cycle processes	EN 62304:2006/AC: 2008 Medical device software - Software life-cycle processes

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Producent:

Guangdong Transtek Medical Electronics Co., Ltd.
Zone A, No.105, Dongli Road, Torch Development District, Zhongshan, 528437, Guangdong, Chiny

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EMC GUIDANCE

- 1) This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.
- 2) Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- 3) **Caution:** This unit has been thoroughly tested and inspected to assure proper performance and operation!
- 4) **Caution:** This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Table 1

Guidance and manufacturer's declaration – electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The device is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	not applicable	
Voltage fluctuations / flicker emissions IEC 61000-3-3	not applicable	

Table 2

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV cont act ± 8 kV air	± 6 kV cont act ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	power supply lines: ± 2 kV input/output lines: ± 1 kV	not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC61000-4-5	line(s) to line(s): ± 1 kV line(s) to earth: ± 2 kV 100 kHz repetition frequency	not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 s	not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50Hz/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

EMC GUIDANCE

Table 3

Guidance and manufacturer's declaration – electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distances: $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.7 GHz where, P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer, d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a) Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device. b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.			

Table 4

Recommended separation distances between portable and mobile RF communications equipment and the device.			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Guangdong Transtek Medical Electronics Co., Ltd.

Address: Zone A, No.105, Dongli Road, Torch Development
District, Zhongshan Guangdong 528437 China

Tel: 86 760 88282980

Fax: 86 760 88282982

Http://www.transtek.cn

E-mail: info@transtek.cn

EU Declaration of Conformity

We Guangdong Transtek Medical Electronics Co., Ltd.

Zone A, No.105, Dongli Road, Torch Development District, 528437

Zhongshan, Guangdong, PEOPLE'S REPUBLIC OF CHINA

declare under our sole responsibility that the product

Product Name: Blood Pressure Monitor

Model name: TMB-1585-BT

Software: V01

Accessories: **Adapter(Model:)KH0601000EW; battery: 4*AAA**

comply with the following directives and regulations:

2014/53/EU (The Radio Equipment Directive)

For the evaluation of the compliance with these Directives and Regulations, the following standards/requirements were applied:

Safety Article3.1(a)	
EMC Article3.1(b)	Draft EN 301 489-1 V2.2.0, Draft EN 301 489-17 V3.2.0
Radio Article3.2	EN 300 328 V2.1.1, Draft EN 303 345 V.1.1.7
Health Article3.1(a)	EN 62479: 2010

The conformity assessment procedure as referenced in Article 17 and detailed in Annex III of the Radio equipment Directive has been followed with the involvement of a notified body:
Notified Body: TUV SUD BABT

Responsible for making this declaration is the:

Manufacturer Authorised representative established within the EU

Person responsible for making this declaration

Print name/Title: Kevin Tan/R&D management

Guangdong Transtek
Medical Electronics Co.,Ltd

2017.9.20

Kevin Tan

SOEHNLE

CONNECT

measure, collect, analyse, advise

Android™: 4.4.4

iOS®: 9.0

Bluetooth®: 4.0

Download on the
App Store

GET IT ON
Google Play

Distributed by
LEIFHEIT AG
Leifheitstraße 1
56377 Nassau/Germany

www.soehnle.com

CE
0123

015358/B 11/2017